

Breaking the Brain's Secret Code



Nathan S. Hicks

and


Larry Sorensen, Advisor

June 07, 2019


The Brain,
The Utah Array,
The Experiment



UTAH ARRAY



A single pyramidal neuron.



Cortex Design


- Thickness: 3.5-mm
- The Utah Array embeds near the middle layer of the cortex.

Electrode Design

- Length: 1.5-mm
- Diameter: 4.0 μm
- Materials: silicone base, platinum tips

Raw Data in the Time Domain:
81 Electrode Voltages vs. Time


Electrode e01 Voltages




Heatmap Plot of All
Electrodes


- 9 x 9 grid
- 81 voltages

A Single Time Sample of Electrode Voltages
vs. Electrode Geometry in Space




Utah Array Electrode Voltages vs. Time





A Single Time Sample of Electrode Voltages vs.
Electrode Geometry in Space




**Heatmap Plot of All
Electrodes**









- 9 x 9 grid
- 81 voltages


Look Up Table (LUT) Transformation Functions



- How to improve the contrast between matrix elements?
- Here, LUT = $f(x)$.



Various Types of Look-Up Table (LUT) Functions and Results Upon Utah Array Electrodes








- [Animation] A time series animation of the Utah Array electrode voltages, in false color.
- [Animation] A time series animation of the Utah Array electrode voltages, in gray scale.


Statistical Analysis:

Histograms, Scatter Plots, and Correlation Coefficients


Electrode e01
Voltages



Probability Density vs. Voltage (mV)
81 Electrodes Using 1×10^6 Time Samples



Electrode e33
Voltages




Voltage Histogram Plots

- Histograms are plotted using all 1,000 seconds of available data.
- Mean of 0.0 for electrodes e01 and e33.


Probability Density Curves

- Mean of 0.0 for all 81 electrode probability density curves.
- A Gaussian analysis is therefore permissible.


Probability Density vs. Voltage (mV)
Time Sample 01




Time Sample 419



Time Sample 976




Histograms for all 81 Electrodes vs. Time




- [Animation] Note the variance in electrode voltages across time.

Scatter Plots

Scatter Plot and Regression Line
Electrode e01 x e01




Scatter Plot and Regression Line
Electrode e01 x e81



Correlation Coefficients, r_{xy}

Correlation Coefficient Histograms
Electrode e01
Bin Size = 1000



Scatter Plot / Correlation


- The variance in X and Y independently.

Correlation Coefficient


- The correlation coefficient, r , is a random sample estimate of both the direction and strength between variables X and Y.
- A positive coefficient result means that as one variable increases, so does the other, and vice versa for a negative coefficient.

$$\begin{aligned}
 r_{xy} &= \frac{\text{cov}(x, y)}{\sqrt{S_x S_y}} = \frac{\sum (x_i - \bar{x})(y_i - \bar{y})}{\sqrt{\sum (x_i - \bar{x})^2 (y_i - \bar{y})^2}} \\
 &= \frac{\text{cov}(x, y)}{\sqrt{\text{var}(y) \text{var}(y)}}.
 \end{aligned}$$

Correlation Coefficient Histograms – Probability Density vs. R-value All Electrodes





- Histogram plots showing the correlation coefficients of each electrode.
- Note that all electrodes have a positive mean and probability density.



- A time series matrix of correlation coefficient heatmap plots.
- Note the generally positive magnitudes of the coefficients.

Correlation Coefficients Across Time




- [Animation] A time series animation of the Utah Array electrode correlation coefficients, in false color.
- [Animation] A time series animation of the Utah Array electrode correlation coefficients, in gray scale.

From Scatter Plots and Correlation Coefficients
To Time Correlation Functions
And Equivalence Classes

Scatter Plot / Covariance Matrix For Two Cases: Electrodes e01–e09 and e28–e36


- The variance in X and Y together.
- 2 type of equivalency classifications are shown between the scatter plot figures.

Electrodes e01–e09




- High correlation in all plots.

Electrodes e28–e36




- Both high and low correlations are observed.
- Increased white space is shown for the highly correlated results.

Scatter Plot / Covariance Matrix Electrodes e01 – e81




Scatter Plot / Covariance Matrix



y-axis Electrodes	x-axis Electrodes Covariance					
	e21-e25	e30-e35	e39-e43	e51-e54	e58-e61	
e14-e17	e21-e25	e30-e35	e39-e43	e51-e54	e58-e61	
e21-e24	e21-e25	e30-e35	e39-e43	e51-e54	e58-e61	
e31-e36	-	e31-e36	-	-	-	
e39-e43	-	e31-e36	e40-e43	-	-	
e49-e54	-	e31-e36	e40-e43	-	-	
e58-e63	-	e31-e36	e40-e43	e50-e53	-	



Equivalency Classifications

- Extracting classification groups via the visible variations in pattern densities.




Cross-Correlation Functions


- Observe the asymmetric characteristic of electrode e33 and those similar in shape.
- Observe the remaining electrode classifications that are symmetric and similar in shape to the auto-correlation.




Auto-Correlation Function of Electrode e01 using 1,000 Seconds of Data




Cross-Correlation Function of Electrodes using 1,000 Seconds of Data Electrodes e01 vs. e33





Auto-Correlation Function of Electrode e33 using 1,000 Seconds of Data



Cross-Correlation Function of Electrodes using 1,000 Seconds of Data Electrodes e33 vs. e01




Explanation of Cross-Correlation Function Characteristics



Sine Wave

'Brain Wave'




<https://youtu.be/L6YJqhbsuFY>

Apply ‘Human Intelligence’ For Pattern Recognition


What Does a Human See?

Lagged Cross-Correlation Peaks - Raw Signals Electrode 1 Lag Frame 230



Human Lag Peak Analysis vs. Electrode Geometry in Space


- How are the lag peaks correlated with the main peak and with each other?




Pattern Recognition – Human Observation

- Determine when 2 or more consecutive frames return nearly identical results.
- Each time step is 1-mSec, making human observable pattern recognition a very challenging method to employ.


Lagged Cross-Correlation Peaks - Raw Signals Electrode 1 Lag Frame 231





Pattern Recognition

- When specific elements of each of the 3 plots give similar cross-correlation magnitudes, the likelihood of electrode signal identification is increased.
- An even stronger likelihood is achieved when the frame-to-frame visualization is nearly identical for the identified electrode signal(s).



Future Work

What Does a Computer See?

Recommendation:

- Apply Artificial Intelligence and Machine Learning algorithm analysis against the Utah Array dataset.
- Machine Learning is particularly good at finding patterns within large datasets.